



MS6110

- **Capacity:** 200 kg / 160 kg
 - **Graduation:** 200 g / 200 g < 100 kg > 500 g
 - **Functions:** BMI , Hold , Tare , Zero
 - **Optional:** ST-3571 / TP2110
- Light weight and easy to carry.
 - Anti slip rubber mat.
 - Flat construction makes stepping on very easy and safe.
 - Easy to measure weight of your baby with 10g accuracy when cuddling your baby.

Weight plays a vital role before, during and after pregnancy.
 Mother's weight during pregnancy directly effects infant's weight.
 Therefore the weight before and during pregnancy and separated weight of mother and infant should be traced carefully and accurately.
 We, Charder, help trace weight of mother and infant via wireless function.

MS6111

- **Capacity:** 250 kg
- **Graduation:** 10 g < 20 kg > 100 g
- **Functions:** BMI , weight measurement , tare , weight lock , dual sided LCD display
- **Optional:** ST-3571 / TP2110



MS 6110




Dual sided LCD display operated by medical professional.



MS 6111

Operated by user him/herself.



ST-3571 Optional carry bag



Features

Rubber mat

Double sided LCD display

Lightweight

MS 6110


Bubble level

Adj. Feet

TP 2110



Flat construction makes stepping on easy and safe.



Weighing up to 250kg.



Measures baby under 20kg by 10g and automatically switch to adult mode above 20kg by 100g.