

# MS7321 Bariatric Floor Scale



## Features

- 250kg (550 lbs) Capacity
- Automatically turns on and weighs when stepped on
- Automatic off
- Reliable, professional-grade measuring
- Wide platform is stable and accommodates bariatric patients
- Low profile allows for easy access
- Self adjusting feet provide prevent tipping
- Non-slip texture



MS7321 Top View

## Specifications

Dimensions	415mm x 360mm x 41 mm (16" x 14" x 1.6" ) WxDxH
Graduation	100g (0.2 lbs)
Power	(4) AA Batteries
Capacity	250 kg (550 lbs)